



The Radiology Information Solution

## Obstetric Ultrasound Module

Standard Forms

**For more details please contact:**

**e:** [sales@hssnet.com](mailto:sales@hssnet.com)

**t:** 01623 489825

HSS, 3rd Floor, i2 Mansfield, Hamilton Court  
Oakham Business Park, Mansfield, NG18 5FB

HSS is part of the  
Wellbeing Software Group



## Copyright

© Healthcare Software Solutions 2014

Registered Office: 3rd Floor i2 Mansfield, Hamilton Court, Oakham Business Park, Mansfield, NG18 5FB

These materials are or comprise restricted and proprietary confidential information of Healthcare Software Solutions. Disclosure to or use by the recipient shall not convey any intellectual property rights in these materials. The right to use these materials by the recipient is subject to restrictions and limitations contained in the Provision of Integrated Care Record System and Associated Services Agreement and related agreements.

## Confidentiality

All information in this document is provided in confidence for the sole purpose of adjudication of the document and shall not be used for any other purpose and shall not be published or disclosed wholly or in part to any other party without HSS prior permission in writing and shall be held in safe custody. These obligations shall not apply to information which is published or becomes known legitimately from some source other than HSS. Many of the product, service and company names referred to in this document are trademarks or registered trademarks. They are all hereby acknowledged.

## Table of Contents

Introduction .....	4
Purpose.....	4
Audience .....	4
Obstetric Reporting Scan Forms.....	5
Early Pregnancy Scan (NW) .....	5
Dating Scan (NW).....	8
Anatomy Scan (NW).....	11
Detailed Scan (NW).....	13
Growth Scan (NW) .....	16
How to Add or Remove a Foetus.....	19
WHEN TO USE 'ADD FOETUS' .....	19
<b>WHEN TO USE 'REMOVE FOETUS'</b> .....	20
Selecting an Alternative Scan Form.....	21
Episodes Screen .....	21
How to Close a Pregnancy or Identify a New Pregnancy .....	22
CLOSE PREGNANCY.....	22
NEW PREGNANCY.....	22
Activating the Obstetric Ultrasound Module.....	23
Obstetric Scan Forms .....	23
Appointment > Exams Per Hospital.....	24
General > Examination Codes .....	24
Obstetric Ultrasound XR Configuration Settings.....	23
Correcting a Pregnancy Record .....	24
CORRECTING OTHER AFFECTED SCANS/ EVENTS .....	25

## Introduction

### Purpose

This CRIB sheet has been written using the default Obstetric Ultrasound Package for the CRIS System specified by the North West Obstetric Clinical Reference Group (CRG).

### Audience

This CRIB sheet is aimed at Obstetricians and Sonographers using the CRIS Obstetric Ultrasound Module, and RIS / PACS System Managers responsible for supporting, configuring and troubleshooting this module.

## Obstetric Reporting Scan Forms

REPORTING AND SONOGRAPHER MODE is recommended for users who are responsible for typing and printing Obstetric reports. CRIS will automatically display / highlight the most appropriate screens and buttons.

When entering details via Obstetric Ultrasound Reports, values can be left blank and the field title will not be printed.

### Early Pregnancy Scan (NW)

To select from the options in each drop down list use [F4] or double click. Select the value you want and [Enter] or double click.

Patient consented to scan

Scan Reason  Method  Uterus

Endometrium  Endometrial thickness  mm

Sac diameter    mm

Gestation Sac  Sac Site  Shape

Yolk Sac  Fetal Pole  Fetal Heartbeat

Adnexa  Free Fluid

#### EARLY PREGNANCY FIELDS (NW)

Patient Consented To Scan - <i>New</i>	A Tick box which allows users to record if the patient has consented to the scan being performed
Scan Reason	Amniocentesis
	New Anomaly
	Booking
	Bleeding
	See Comments
	CVS
	Diabetic
	Dating & NT - <i>New</i>
	Dating
	Fetal Study
	Growth
	Medical Disorder
	NT - New
	Prev Fetal Abnorm
	Poor Obs History
Pain	
Prev Obs History	
Reassurance	


	Rescan Anomaly
	Transfer Booker
	Twins
	Viability
Method	TA + TV
	Trans Abdominal
	Trans Vaginal
Uterus	Bulky
	Normal
Endometrium	Abnormal
	Normal
	Thickened
	Thin
Endometrial thickness (mm)	Enter in mm.
Sac diameter <input type="text" value="0.0"/> <input type="text" value="0.0"/> <input type="text" value="0.0"/> mm	Enter measurements in mm. in all fields.
Gestation Sac	More Than One
	Not Present
	Present
Sac Site	Abnormal
	Cervical
	Cornual
	Ectopic
	Normal
Shape	Irregular
	Regular
Yolk Sac	See Comments / Seen / Not Seen
Fetal Pole	
Fetal Heartbeat	More Than One
	Not Present
	Present
Adnexa	See Comments / Seen / Not Seen
Free Fluid	

### EARLY PREGNANCY MEASUREMENTS

**CRL**  mm  Wks  Days 

**Gestational Age By Scan**  Wks  Days **EDD**

**EARLY PREGNANCY MEASUREMENTS WITH LMP (IF ENABLED)**

CRL	<input type="text" value="9.5"/>	mm	<input type="text" value="7"/>	Wks	<input type="text" value="0"/>	Days	
LMP	<input type="text" value="25/12/2012"/>		<input type="text" value="7"/>	Wks	<input type="text" value="0"/>	Days	EDD <input type="text" value="01/10/2013"/>
Gestational Age By CRL			<input type="text" value="7"/>	Wks	<input type="text" value="0"/>	Days	EDD <input type="text" value="01/10/2013"/>

Enter CRL (mm) and press [Enter] to automatically calculate the Gestational Age by CRL (Crown Rump Length) and EDD (Estimated Date of Delivery). These calculations are based on Robinson key equations and calculations.

The Gestational Age by Scan ,or EDD can be manually overridden by the User but this will be displayed as follows:

Gestational Age By USER  Wks  Days EDD

Please note: If the EDD date is removed at this point the EDD on future scan forms will have to be manually entered. Similarly if the EDD is manually amended at this point the EDD, future scans will also display this amended date as it is assumed that the 'User' has chosen to override the automatic calculation. The only exception to these rules is if the Trust has enabled the EDDCutOffDays XR Setting which defines the no.days before an EDD is considered accurate (i.e. 56 Days – 8 weeks).

Screening Incomplete? <input type="checkbox"/>	<b>New</b> - Incomplete screenings can be recorded by clicking the 'Screening Incomplete?' tick box.
Comments:	The comments section should be used to enter additional information / details regarding the scan. To type additional comments left click beneath <b>Comments:</b> and begin typing.

Please note: Obstetric forms can also be configured not to display LMP, print EDD or suggested rescan date.

## Dating Scan (NW)

Patient consented to scan

Scan Reason  Scan Method  Fetal Heartbeat   
 Adnexa  Placenta  Chorionicity   
 NT  mm Sac Diameter

### DATING SCAN (NW)

Patient Consented to Scan	Tick box to record if the patient has consented to the scan being performed.
Scan Reason	Amniocentesis
	New Anomaly
	Booking
	Bleeding
	See Comments
	CVS
	Diabetic
	Dating & NT - <i>New</i>
	Dating
	Fetal Study
	Growth
	Medical Disorder
	NT - New
	Prev Fetal Abnorm
	Poor Obs History
	Pain
	Prev Obs History
Reassurance	
Rescan Anomaly	
Transfer Booker	
Scan Method	TA + TV
	Trans Abdominal
	Trans Vaginal
Fetal Heartbeat	More Than One
	Not Present
	Present
Adnexa	See Comments/ Seen/ Not Seen
Placenta	Anterior
	Fundal
	Left Lateral
	Posterior
	Right Lateral

Chronicity	MA / DC
	DA / DC
	Monoamniotic / Dichorionic – <i>New</i>
	MA / MC
NT (mm)	Nuchal Translucency entered in mm
Sac Diameter (mm)	Enter in mm.

#### DATING MEASUREMENTS

CRL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD	<input type="text"/>
Anatomy	<input type="text" value="20"/> Wks	<input type="text" value="0"/> Days	Date	<input type="text" value="03/08/2013"/>
Suggested Rescan Date	<input type="text" value="34"/> Wks	<input type="text" value="0"/> Days	Date	<input type="text" value="09/11/2013"/>

#### DATING MEASUREMENTS WITH LMP (IF ENABLED)

CRL	<input type="text" value="56.0"/> mm	<input type="text" value="12"/> Wks	<input type="text" value="1"/> Days	
HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
LMP	<input type="text" value="25/12/2012"/>	<input type="text" value="11"/> Wks	<input type="text" value="4"/> Days	EDD <input type="text" value="01/10/2013"/>
Gestational Age By CRL	<input type="text" value="11"/> Wks	<input type="text" value="4"/> Days	EDD	<input type="text" value="01/10/2013"/>
Anatomy	<input type="text" value="20"/> Wks	<input type="text" value="0"/> Days	Date	<input type="text" value="14/05/2013"/>
Suggested Rescan Date	<input type="text" value="34"/> Wks	<input type="text" value="0"/> Days	Date	<input type="text" value="20/08/2013"/>

Enter CRL, HC and FL as applicable and press [Enter] to automatically calculate the Gestational Age by CRL (Crown Rump Length), or HC if this is not entered or FL if CRL and HC have not been entered.

For the Dating Scan the EDD is carried from the previous scan as applicable, if this is not present then EDD calculated on CRL, if this is not present then EDD is calculated on HC based on Robinson and Chitty key equations / calculations respectively.

The Gestational Age by Scan ,or EDD can be manually overridden by the User but this will be displayed as follows:

Gestational Age By USER  Wks  Days EDD

The Gestational Age by Scan can also be overridden to reflect a particular measurement – i.e HC or CRL by matching the Gestational Age by scan Wks and Days to the appropriate measurement.

Gestational Age By HC  Wks  Days EDD

Please note: If the EDD date is removed at this point the EDD on future scan forms will have to be manually entered. Similarly if the EDD is manually amended at this point the EDD, future scans will also display this amended date as it is assumed that the 'User' has chosen to override the automatic calculation.

Screening Incomplete?

**New** - Incomplete screenings can be recorded by clicking the 'Screening Incomplete?' tick box.

Comments:

The comments section should be used to enter additional information / details regarding the scan. To type additional comments left click beneath **Comments:** and begin typing.

Please note: The Anatomy date field is only designed as a guide to assist in planning future appointments. This defaults to 20 weeks but can be changed to allow you to calculate any date. It is not stored on the report and is only accurate once the EDD has been entered. Forms can also be configured not to display LMP, print EDD or suggested rescan date.

## Anatomy Scan (NW)

Patient consented to scan <input type="checkbox"/>	Scan Reason <input type="text"/>	
Fetal Heartbeat <input type="text"/>	4C Heart <input type="text"/>	Outflow <input type="text"/>
Cranial Structures <input type="text"/>	Cerebellum <input type="text"/>	Face <input type="text"/>
Stomach/Diaphragm <input type="text"/>	Kidneys <input type="text"/>	Bladder <input type="text"/>
Abdo Wall <input type="text"/>	Spine <input type="text"/>	
Upper Limbs <input type="text"/>	Lower Limbs <input type="text"/>	
Placenta <input type="text"/>	Site <input type="text"/>	Liquor Volume <input type="text"/>

ANATOMY SCAN	
Patient Consented to scan	Tick box to record if the patient has consented to the scan being performed.
Fetal Heartbeat	More Than One
	Not Present
	Present
4C Heart	See Comments/ Seen/ Not Seen
Cranial Structures	
Cerebellum	
Face	
Stomach/ Diaphragm	
Kidneys	
Bladder	
Abdo Walls	
Spine	
Upper Limbs	
Lower Limbs	
Placenta	Anterior
	Fundal
	Left Lateral
	Posterior
Site	Abutting OS
	Clear of OS
	Covering OS
	Low, clear of OS
	Low Lying
	Overlying OS
	See Comments
Touching OS	

Liquor Volume (cm)	Low – Additional fields will appear to complete AFI values. AFI <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> = <input type="text" value="0.0"/>
	Normal
	Oligohydramnios – Additional fields will appear to complete AFI values.
	Polyhydramnios – Additional fields will appear to complete AFI values.
	Raised- complete – Additional fields will appear to complete AFI values.

### ANATOMY MEASUREMENTS

HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
AC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD <input type="text"/>	

### ANATOMY MEASUREMENTS WITH LMP (IF ENABLED)

HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
AC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
LMP	<input type="text" value="25/12/2012"/>	<input type="text" value="19"/> Wks	<input type="text" value="4"/> Days	EDD <input type="text" value="01/10/2013"/>
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD <input type="text"/>	

Enter HC, AC, and FL as applicable and press [Enter] to automatically calculate the Gestational Age by HC (Head Circumference) or FL if this is not entered.

For the Anatomy Scan the EDD is carried from the previous scan, if this is not present then EDD calculated on HC, if this is not present then EDD calculated on FL. These are based on Chitty key equations / calculations respectively.

The Gestational Age by Scan ,or EDD can be manually overridden by the User but this will be displayed as follows:

Gestational Age By USER  Wks  Days EDD

The Gestational Age by Scan can also be overridden to reflect a particular measurement – i.e HC or FL by matching the Gestational Age by scan Wks and Days to the appropriate measurement.

Gestational Age By HC  Wks  Days EDD

Please note: If the EDD date is removed at this point the EDD on future scan forms will have to be manually entered. Similarly if the EDD is manually amended at this point the EDD, future scans will also display this amended date as it is assumed that the ‘User’ has chosen to override the automatic calculation.

Screening Incomplete? <input type="checkbox"/>	<b>New</b> - Incomplete screenings can be recorded by clicking the ‘Screening Incomplete?’ tick box.
Comments:	The comments section should be used to enter additional information / details

regarding the scan. To type additional comments left click beneath **Comments:** and begin typing.

### Detailed Scan (NW)

Patient consented to scan  Scan Reason

Fetal Heartbeat <input type="text"/>	4C Heart <input type="text"/>	Outflow <input type="text"/>
Ventricles <input type="text"/>	CSP <input type="text"/>	Cerebellum <input type="text"/>
Face <input type="text"/>	Lips <input type="text"/>	Profile <input type="text"/>
Stomach/Diaphragm <input type="text"/>	Kidneys <input type="text"/>	Bladder <input type="text"/>
Spine Sagittal <input type="text"/>	Coronal <input type="text"/>	Transverse <input type="text"/>
Upper Limbs <input type="text"/>	Hands <input type="text"/>	Abdo Wall <input type="text"/>
Lower Limbs <input type="text"/>	Feet <input type="text"/>	Cord <input type="text"/>
Placenta <input type="text"/>	Site <input type="text"/>	Liquor Volume <input type="text"/>
Fetal Movement <input type="text"/>		
EDF <input type="text"/>	PI <input type="text" value="0.0"/>	RI <input type="text" value="0.0"/>
		SD <input type="text" value="0.0"/>

DETAILED SCAN (NW)	
Patient Consented to scan	Tick box to record if the patient consented to the scan being performed.
Fetal Heartbeat	More Than One
	Not Present
	Present
4C Heart	See Comments / Seen / Not Seen
Outflow	
Ventricles	
CSP (Cavum Septum Pellucidum)	
Cerebellum	
Face	
Lips	
Profile	
Stomach/ Diaphragm	
Kidneys	
Bladder	
Spine Sagittal	
Coronal	
Transverse	
Upper Limbs	
Hands	
Abdo Wall	

Lower Limbs	
Feet	
Cord	2 Vessels
	3 Vessels
	Not Seen
Placenta	Anterior
	Fundal
	Left Lateral
	Posterior
	Right Lateral
Site	Abutting OS
	Clear of OS
	Covering OS
	Low, clear of OS
	Low Lying
	Overlying OS
	See Comments
	Touching OS
Liquor Volume (cm)	Low – Additional fields will appear to complete AFI values. AFI <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> = <input type="text" value="0.0"/>
	Normal
	Oligohydramnios – Additional fields will appear to complete AFI values.
	Polyhydramnios – Additional fields will appear to complete AFI values.
	Raised- complete – Additional fields will appear to complete AFI values.
Fetal Movement	See Comments/ Seen/ Not Seen
EDF - End Diastolic Flow	Absent / Present / Reversed
PI - Pulsatility Index	Enter as required
RI - Resistivity Index	
S/D - Systolic / Diastolic Ratio	

**DETAILED MEASUREMENTS**

HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
AC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
Cerebellum	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD <input type="text"/>	

**DETAILED MEASUREMENTS WITH LMP (IF ENABLED)**

HC	<input type="text" value="0.0"/>	mm	<input type="text"/>	Wks	<input type="text"/>	Days	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
AC	<input type="text" value="0.0"/>	mm	<input type="text"/>	Wks	<input type="text"/>	Days	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
FL	<input type="text" value="0.0"/>	mm	<input type="text"/>	Wks	<input type="text"/>	Days	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Cerebellum	<input type="text" value="0.0"/>	mm	<input type="text"/>	Wks	<input type="text"/>	Days	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
LMP	<input type="text" value="25/12/2012"/>		<input type="text" value="19"/>	Wks	<input type="text" value="4"/>	Days	EDD	<input type="text" value="01/10/2013"/>		
Gestational Age By Scan	<input type="text" value="0"/>		<input type="text" value="0"/>	Wks	<input type="text"/>	Days	EDD	<input type="text"/>		

Enter HC, AC, FL and Cerebellum - **New** as applicable and press [Enter] to automatically calculate the Gestational Age by HC (Head Circumference) or FL if this is not entered.

For the Detailed Scan the EDD is carried from the previous scan, if this is not present then EDD calculated on HC, if this is not present then EDD calculated on FL. These are based on Chitty key equations / calculations respectively.

The Gestational Age by Scan ,or EDD can be manually overridden by the User but this will be displayed as follows:

Gestational Age By USER  Wks  Days EDD

The Gestational Age by Scan can also be overridden to reflect a particular measurement – i.e FL by matching the Gestational Age by scan Wks and Days to the appropriate measurement.

Gestational Age By FL  Wks  Days EDD

**Please note:** If the EDD date is removed at this point the EDD on future scan forms will have to be manually entered. Similarly if the EDD is manually amended at this point the EDD, future scans will also display this amended date as it is assumed that the ‘User’ has chosen to override the automatic calculation.

Screening Incomplete? <input type="checkbox"/>	<b>New</b> - Incomplete screenings can be recorded by clicking the ‘Screening Incomplete?’ tick box.
Comments:	The comments section should be used to enter additional information / details regarding the scan. To type additional comments left click beneath <b>Comments:</b> and begin typing.

## Growth Scan (NW)


Patient consented to scan <input type="checkbox"/>	Scan Reason <input type="text"/>		
Fetal Heartbeat <input type="text"/>	Presentation <input type="text"/>	Fetal Movements <input type="text"/>	
Stomach / Diaphragm <input type="text"/>	Kidneys <input type="text"/>	Bladder <input type="text"/>	
Placenta <input type="text"/>	Site <input type="text"/>	Liquor Volume <input type="text"/>	
EDF <input type="text"/>	PI <input type="text" value="0.0"/>	RI <input type="text" value="0.0"/>	S/D <input type="text" value="0.0"/>

### GROWTH SCAN (NW)


Patient Consented to scan	Tick box to record if the patient consented to the scan being performed.
Fetal Heartbeat	More Than One
	Not Present
	Present
Presentation	Breech
	Breech Legs Extended
	Footling Breech
	Breech Maternal Left
	Breech Maternal Right
	Cephalic
	Cephalic Maternal Left
	Cephalic Maternal Right
	Head to Mat. Left
	Oblique Lie
	Heat to Mat. Right
	Transverse Lie
	Vertical
Fetal Movements	See Comments/ Seen/ Not Seen
Stomach/ Diaphragm	
Kidneys	
Bladder	
Placenta	Anterior
	Fundal
	Left Lateral
	Posterior
	Right Lateral
Site	Abutting OS
	Clear of OS
	Covering OS
	Low, clear of OS
	Low Lying
	Overlying OS

	See Comments
	Touching OS
Liquor Volume (cm)	Low – Additional fields will appear to complete AFI values. AFI <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> + <input type="text" value="0.0"/> = <input type="text" value="0.0"/>
	Normal
	Oligohydramnios – Additional fields will appear to complete AFI values.
	Polyhydramnios – Additional fields will appear to complete AFI values.
	Raised- complete – Additional fields will appear to complete AFI values.
EDF - End Diastolic Flow	Absent / Present / Reversed
PI - Pulsatility Index	Enter as required
RI - Resistivity Index	
S/D - Systolic / Diastolic Ratio	

### GROWTH MEASUREMENTS

HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
AC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD <input type="text"/>	
EFW by	<input type="text"/>	<input type="text" value="0.0"/> gms		

### GROWTH MEASUREMENTS WITH LMP (IF ENABLED)

HC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
AC	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
FL	<input type="text" value="0.0"/> mm	<input type="text"/> Wks	<input type="text"/> Days	
LMP	<input type="text" value="25/12/2012"/>	<input type="text" value="31"/> Wks	<input type="text" value="4"/> Days	EDD <input type="text" value="01/10/2013"/>
Gestational Age By Scan	<input type="text" value="0"/> Wks	<input type="text" value="0"/> Days	EDD <input type="text"/>	
EFW by	<input type="text"/>	<input type="text" value="0.0"/> gms		

Enter HC, AC, and FL as applicable and press [Enter] to automatically calculate the Gestational Age by HC (Head Circumference) or FL if this is not entered.

For the Growth Scan the EDD is carried from the previous scan, if this is not present then EDD calculated on HC, if this is not present then EDD calculated on FL. These are based on Chitty key equations / calculations respectively.

The Gestational Age by Scan ,or EDD can be manually overridden by the User but this will be displayed as follows:

Gestational Age By USER  Wks  Days EDD

The Gestational Age by Scan can also be overridden to reflect a particular measurement – i.e FL by matching the

Gestational Age by scan Wks and Days to the appropriate measurement.

Gestational Age By FL  Wks  Days EDD

Please note: If the EDD date is removed at this point the EDD on future scan forms will have to be manually entered. Similarly if the EDD is manually amended at this point the EDD, future scans will also display this amended date as it is assumed that the 'User' has chosen to override the automatic calculation.

EFW - Estimated Fetal Weight	EFW by <input type="text"/>	<input type="text" value="0.0"/> gms
	C	Campbell (AC only)
	H1	Hadlock AC + FL
	H2	Hadlock HC + AC + FL
	S	Shepherd AC + BPD ( <i>Legacy field value for previous data integrity</i> )
<p>This <input type="text" value="0.0"/> gms field will be completed automatically based on the combination of AC and FL measurements, or according to EFW priority specified via the XR Setting. This can also be overridden via the 'EFW by' field by pressing [F4] and selecting from the resulting list.</p> <p>The default configuration of CRIS is to calculate the EFW in grams using the equation priority of Hadlock 1 - If the AC and FL are entered or Campbell if only AC is entered. It is however possible to prevent AC only calculations by setting the XR setting REPORT.ObsAllowEFWbyACOnly = N so that this equation is ignored.</p> <p>CRIS can also be configured using the XR Setting REPORT.ObsEFWPriority to change the order of equation priority from the default. Changing the default priority order will also allow the users to choose from the list of options using [F4] on the EFW field to calculate EFW on the report. <u>Please note:</u> Unless otherwise stated HSS recommend H2,H1 (Hadlock 2 then Hadlock 1).</p>		
Screening Incomplete? <input type="checkbox"/>	<b>New</b> - Incomplete screenings can be recorded by clicking the 'Screening Incomplete?' tick box.	
Comments:	The comments section should be used to enter additional information / details regarding the scan. To type additional comments left click beneath <b>Comments:</b> and begin typing.	

## Additional Functionality

### How to Add or Remove a Foetus

REPORT E-21914 Attended-03/04/2008

Fetus/Sac 1 of 1 Early Scan (North West) Start 03/03/2008 End Unknown Key 21914

Scan Reason  Method  Uterus

Endometrium  Endometrial thickness

Sac diameter    mm

Gestation Sac  Sac Site

Yolk Sac  Fetal Pole

Adnexa  Free Fluid

CRL  mm  Wks

Gestational Age By Scan  Wks


Comments:


- Undo
- Redo
- Verify
- Mark as unchecked (Check)
- Add addendum
- Set Actioned by AMTRAIN
- Load PACS images for Exam
- Add foetus**
- Remove foetus Add a foetus
- Close Pregnancy
- New Pregnancy
- Change Scan Form
- Edit Series
- Properties
- Restricted

### WHEN TO USE 'ADD FOETUS'

To report on more than one fetus (e.g. twins) you should right-click via the Report Editor and select 'Add foetus' from the resulting menu. This will load another scan form into the same report and display 'Fetus/Sac 1 of 2' on the current, and all future scan forms. For triplets/ quadruplets etc. continue to select 'Add foetus' until the correct number of scan forms have appeared.

Double click a graph to enlarge the graph and double click an enlarged graph to return back to normal view with all graphs.






Each Fetus is displayed using a different colour (shown via the key at the top) and shape to identify the different measurements for each fetus. The colour and shape is displayed on screen, and printed as unique shapes via black and white printouts.

#### WHEN TO USE 'REMOVE FOETUS'

'Remove Foetus' is designed to identify that a 'FETUS IS NO LONGER VIABLE' and can be selected via the right-click menu via the Report Editor at the point you notice this is the case. The Fetus details will remain on the current report but will not appear on any future scans. This tool is not designed to close an event with only one fetus, if there is only one fetus please use 'Close Pregnancy'.

**Please note:** If you mistakenly select 'Add Foetus' and have not saved the report you should select [F5] to clear the report and the incorrectly added fetus. If the report has already been saved you will have to use 'Remove Foetus' - the 2<sup>nd</sup> fetus will remain on the original report, but not on any subsequent scans.

## Selecting an Alternative Scan Form



Selecting 'Change Scan Form' allows the user to change scan forms without going back to the Event Details screen and changing the exam code – i.e. changing a 'Detailed Scan' form to an 'Anomaly Scan' form. However it is essential that you select 'Change Scan form' and load the new scan form BEFORE any details are entered in order to maintain data integrity.

## Episodes Screen

All scans related to each pregnancy term are visible via the Episodes screen tab. The folder for each pregnancy term can be opened to display all associated scans, and the number of Fetus and EDD are displayed on the folder.



## How to Close a Pregnancy or Identify a New Pregnancy

REPORT E-19214 Attended-12/03/2007

Fetus/Sac 1 of 1 Early Scan (North West) Start 21/01/2008 End 28/10/2007 Key 19214

Scan Reason

Endometrium  Endometrial thick

Sac diameter    mm

Gestation Sac

Yolk Sac

Adnexa

CRL  mm

Gestational Age By CRL

Method

Sac

Fetal

Free

Uterus

Shape

Heartbeat

- Undo
- Redo
- Verify
- Mark as unchecked (Check)
- Add addendum
- Set Actioned by AMTRAIN
- Load PACS images for Exam
- Add foetus
- Remove foetus
- Close Pregnancy
- New Pregnancy**
- Change Scan Form New Pregnancy
- Edit Series
- Properties
- Restricted
- Delete Report
- Print

Comments:

### CLOSE PREGNANCY

‘Close Pregnancy’ is used to identify that a pregnancy is not going to go to full term. Consequently, this should only be used in the event of a scan which results in identifying the fetus is no longer viable, or an anomaly scan which will end with a termination.

If the next scan takes place at least one month after using ‘Close Pregnancy’ CRIS will automatically allocate a new pregnancy key and EDD – otherwise it will be necessary to use ‘New Pregnancy’ and manually correct the EDD.

Please note: ‘Close Pregnancy’ should NOT be used at the time of booking / reporting a new pregnancy. **NEW** - A new security setting of REPORT.CLOSE\_PREGNANCY has been created in order to control whether the ‘close pregnancy’ option is available or not.

### NEW PREGNANCY

If a previous pregnancy does not go to full term (i.e. in the event of miscarriage or termination), and the patient becomes pregnant again during the original/previous pregnancy term, it will be necessary to use ‘New Pregnancy’. This should only be used when booking/reporting the first scan of the new pregnancy term.

It will however be necessary to check / manually correct the EDD of this new pregnancy.

Please note: If ‘New Pregnancy’ is not used, CRIS will assume that the scan is part of previous pregnancy term and use the same EDD, pregnancy key and Episode.

## Obstetric Ultrasound Module Access Settings

### ACTIVATING THE OBSTETRIC ULTRASOUND MODULE

To begin using the CRIS Obstetric Ultrasound Module the RIS System Administrators will need to assign the appropriate Obstetric forms to an examination code using the 'Scan Reason' field located in the CRIS System Tables via either the **EXAMINATION TABLE** (this is SHA / Consortium wide) or **EXAMS PER HOSPITAL TABLE** which is site specific. It will also be necessary to review/configure the relevant XR Settings in context to Trust working practice.

### SECURITY SETTINGS

There is only one security settings associated with the Obstetric Ultrasound Module and that is **REPORT.CLOSE\_PREGNANCY** which controls access to the setting via the right-click menu in the report editor and will need to be assigned to all applicable user groups (roles).

## Obstetric Module Table Configuration

### OBSTETRIC ULTRASOUND XR CONFIGURATION SETTINGS

Please refer to CRIS\_CRIB\_IM\_354\_Obstetric\_Ultrasound\_Configuration.doc for full details of all configuration settings which enable the Trust to customise the module as applicable.

### OBSTETRIC SCAN FORMS

The following scan forms are available via the *Standard Obstetric Ultrasound Forms* and can be assigned to appear automatically when reporting a specific Obstetric examination. This configuration is undertaken via either Exams Per Hospital (Site level configuration) or Examination codes (System Wide configuration).

SCAN FORM CODE	NAME/DESCRIPTION
ENW	Early pregnancy Scan (1 <sup>st</sup> Trimester Viability Scan)
DATNW	Dating Scan (1 <sup>st</sup> Trimester Scan)
ANATNW	Anatomy Scan (Basic 2 <sup>nd</sup> Trimester Scan)
DETNW	Detailed Scan (An alternative and more detailed 2 <sup>nd</sup> Trimester Scan)
GROWNW	Growth Scan (Late 2 <sup>nd</sup> and general 3 <sup>rd</sup> Trimester Scan)
FOR USE WITH FETAL ANOMALY MODULE IF PURCHASED	
ANOMNW	Fetal Anomaly Scan Form
BIOPNW	Biophysical Profile Scan Form
OUTCNW	Pregnancy Outcome Form (Post Natal)

## APPOINTMENT > EXAMS PER HOSPITAL

HSS recommend that Obstetric 'Scan Reasons' are primarily applied via **EXAMS PER HOSPITAL** as 'Forms' assigned at this level are 'Site' specific rather than SHA / Consortium to enable each Trust to make their own decision to which forms should appear when using different Examination codes. Each Trust should however ensure a consistent approach across all sites to negate any compatibility/data integrity issues between sites (i.e. if an early pregnancy was performed at Hospital A, but the Dating Scan at Hospital B).

## GENERAL > EXAMINATION CODES

Setting Scan Forms via GENERAL > EXAMINATION CODES is typically only available to Locally Deployed CRIS Systems not Datacentre models as 'Scan Reasons' assigned at this level will be applicable to the entire system / Consortium. Datacentre / Consortium customers should not complete 'Scan Reasons' here unless specifically advised to do so by the Consortium focus group, or CRIS helpdesk.

# Troubleshooting

## CORRECTING A PREGNANCY RECORD

A typical indication that the user has not selected [New Pregnancy] would be if the EDD has been carried over from the previous pregnancy term.

To confirm this go to the [Episodes] folder and in normal circumstances there should be a new pregnancy folder with a new pregnancy key and the correct EDD. However, if this is not visible it is likely that the most recent scan has been incorrectly attached to a previous pregnancy term. Consequently by opening the last pregnancy folder and reviewing the contents you should immediately notice the most recent scan is actually in this folder.

To correct this and create a new Episode 'Pregnancy Folder' and EDD - Load the report for the first scan (Event) of the new pregnancy term


Right click and choose [New Pregnancy] followed by [Finished] and [Save] to allocate a new pregnancy key and amend EDD to the correct date for the new pregnancy term.

Reload the report and confirm that the Pregnancy key (which is displayed on the top right of the report) has changed to match the Event key (which is in blue text on the top left of the screen).

You should now check there is a new pregnancy folder with the same pregnancy key and EDD displayed via the Episodes screen.

**CORRECTING OTHER AFFECTED SCANS/ EVENTS**

If there are any subsequent scans for this pregnancy you will also need to manually allocate the same new pregnancy key and EDD. Consequently you should go to Events screen and right click the next affected scan and chose [Properties].



Replace the incorrect [Preg Key] with the new [Preg Key] / Event key followed by [Save].

Once the Pregnancy Key is revised you will need to reload the report and check/amend the measurements or EDD. To correct Obstetric Measurements or EDD in a report it may first be necessary to [Unverify] the original report, edit the relevant details and re-verify the report.

It is then then strongly recommend that an Addendum should be entered to stats "Please be advised that the measurements and EDD have been updated following the initial issue of this report" or words to that effect before re-verifying the report.

This process will need to be undertaken by a Key Users (i.e. Superintendent Sonographers or RIS / PACS Manager) as Unverify permission should not be given to all users as it would represent a clinical risk to unverify / edit / re-verify if this process is not then supported by an addendum entry.

Repeat this process for any other scans/ Events.

**Glossary of Terms**

AC	Abdominal Circumference
AFI	Amniotic Fluid Index
CRL	Crown Rump Length
EDF	End Diastolic Flow
FL	Femur Length
EDD	Estimated Delivery Date
EFW	Estimated Fetal Weight
HC	Head Circumference
LMP	Last Menstrual Period
NC	Nuchal Translucency

PI	Pulsatility Index
RI	Resistivity Index
S/D	Systolic / Diastolic Ratio

## Associated Documentation

**CRIS\_CRIB\_IM\_354\_Obstetric Ultrasound Configuration.doc** – Details of all XR Settings which enable the Trust to customise the module as applicable.

**CRIS\_PD118\_CRIS\_Obstetrics\_Module\_Calcs\_NW.doc** – Technical Document citing the key equations and calculation logic for the Obstetrics module.

## Document Control

<b>Title</b>	Obstetric Ultrasound Reporting – Standard Forms (North West)		
<b>Owner</b>	HSS Training Manager	<b>Date Created</b>	01/09/2009
<b>File Ref.</b>	CRIS_CRIB_IM_354_Obstetric_Reporting_V2.1.doc		
<b>CRIS Version</b>	2.09.10p		
<b>Change History</b>			
<b>Issue</b>	<b>Date</b>	<b>Author / Editor</b>	<b>Details of Change</b>
1.0	01/09/2009	David Costin	First Issue
1.1	04/07/2011	David Costin	Amendments to version number for 2.09.10e
1.2	04/11/2011	David Costin	Amendments to facilitate changes made in 2.09.10i
D2.0	30/09/2013	Emma Savage-Mady	Amendments to facilitate changes made from 2.09.10k to 2.09.10p7.
V2.0	01/06/2014	Emma Savage-Mady	Version V2.0 for release 2.09.10p including name change to Healthcare Software Solutions.
V2.1	05/01/2017	Emma Savage-Mady	Amendments to underline use of Addendum in correcting measurements.
<b>Review Date</b>	01/09/2015		