



## Core CRIS Augnito Voice Commands

All voice commands for use in conjunction with Core CRIS reporting.

| Microphone                                   |  |
|--|--|
| Microphone on                                | Click the Mic Icon on screen or use Speech Mike Record Button.<br><b>Please note:</b> This command cannot be implemented with Augnito as there is no concept of sleep mode / background listening. |
| Microphone off<br>Mic stop<br>Dictation stop | Switches microphone off.   |
| Wake up                                      | Click the Mic Icon on screen or use Speech Mike Record Button.<br><b>Please note:</b> This command cannot be implemented with Augnito as there is no concept of sleep mode / background listening. |
| Go to sleep                                  | Click the Mic Icon on screen or use Speech Mike Record Button.<br><b>Please note:</b> This command cannot be implemented with Augnito as there is no concept of sleep mode / background listening. |
| Insert                                       |  |
| Insert title                                 | Insert a title in the current section.   |
| Insert all titles                            | Insert clinical history and all exam titles, and press Enter for each.   |
| Insert exam titles                           | Insert all exam titles, and press Enter for each.  |
| Start report                                 | Insert exam title in the current section and moves cursor to the next line.  |
| Next exam start report                       | Combines the 'next exam' and 'start report' commands.  |
| Punctuate                                    |  |
| Comma  | Add a comma.   |
| Full stop                                    | Add a full stop.   |
| Open bracket<br>Close bracket                | Add a left / right bracket (parenthesis).  |
| Save   |  |
| Save report                                  | Save a report as unverified.   |
| Finish report                                | Press the Finished button to go to the Finished screen, where you can mark for checking.   |

|  |  |
|--|--|
| Clear patient                              | Close the patient record without saving (press F5).                                  |
| Verify report                              | Save a report as verified.   |
| <b>Show</b>                                |  |
| PACS Off                                   | Turns off PACS auto image loading.   |
| PACS On                                    | Turns on PACS auto image loading.  |
| Load Image                                 | Sends PACS message to load images for this exam.                                     |
| Show request card<br>Hide request card     | Show / hide a scanned request card in the report editor.                             |
| <b>Navigate</b>                            |  |
| New line                                   | Move the cursor down one line.   |
| New paragraph<br>Paragraph<br>New para     | Start a new paragraph by moving the cursor down two lines.                           |
| First field                                | Jump to the first [field] in a report, indicated by [square brackets].               |
| Next field                                 | Jump to the next [field] in a report, indicated by [square brackets].                |
| Previous field                             | Jump to the previous sequential [field] in a report, indicated by [square brackets]. |
| Last field                                 | Jump to the final [field] in a report, indicated by [square brackets].               |
| Next exam                                  | Jump to next exam section.   |
| Previous exam                              | Jump to previous exam section (or Clinical History section).                         |
| Got to top<br>Go to bottom                 | Move cursor to bottom or top of the current exam field.                              |
| Go to end of line<br>Go to line end        | Move cursor to the end of this line.   |
| Go to start of line<br>Go to line start    | Move cursor to the start of this line.   |
| Go to document start<br>Go to document end | Move cursor to the start of the report.<br>Move cursor to the end of the report.     |
| Move up / down <n> lines                   | Move cursor up / down.   |

| Select                       |  |
|------------------------------|--|
| Select previous word         | Selects the word before the current cursor position.<br><b>Please note:</b> This only works if there is no space between the last word and current cursor position.    |
| Select [word]                | Searches backwards from the current cursor position and selects the first instance of the [word] found.  |
| Select [word] through [word] | Select a text range. This phrase must be used even if selecting only two words. e.g. "Select clinical through history" and not "select clinical history"               |
| Select again                 | Select the previous instance of a selected word. Repeat to select other instances of that word.  |
| Select unrecognised word     | Select unrecognised text (indicated by underlined space).  |
| Select next [0 - 9] + words  | Select the next [0 - 9] + words.   |
| Unselect that                | Deselect selected text.  |
| Edit                         |  |
| Cut it                       | Cut selected text and save to the clipboard.   |
| Copy it                      | Copy selected text and save to the clipboard.  |
| Paste it                     | Paste clipboard text.  |
| Undo it<br>Undo that         | Undo last action (for example unintentional deletion).   |
| Delete                       |  |
| Delete that                  | Delete selected text.  |
| Scratch that                 | Delete all text that was dictated since your last pause.   |
| Delete line                  | Delete the current line.   |
| Delete last <n> words        | Delete last number of words.   |
| Delete previous word         | Delete the word prior to the current cursor location.  |
| Format                       |  |
| All caps on / off            | All text upper-case until you turn off.  |
| Capitalise that              | Capitalise selected text.  |
| Capitalise last word         | Capitalise the word before the current cursor position.<br><b>Please note:</b> This only works if there is no space between the last word and current cursor position. |

|                                  |   |
|----------------------------------|---|
| Begin bold, Start bold           | Start bold text.  |
| End bold, Stop bold              | Stop bold text.   |
| Bold it, Bold that               | Bold selected text.   |
| Begin italic/s, Start italic/s   | Start italic text.  |
| End italic/s, Stop italic/s      | Stop italic text.   |
| Italicise it, Italicise that     | Italicise selected text.  |
| Begin underline, Start underline | Start underlining text.   |
| End underline, Stop underline    | Stop underlining text.  |
| Underline it, Underline that     | Underline selected text.  |
| Underline last word              | Underline the word before the current cursor position.<br><b>Please note:</b> This only works if there is no space between the last word and current cursor position. |

## Document Control

|                       |  |                   |  |
|-----------------------|--|-------------------|--|
| Title                 | Augnito Voice Commands                           |                   |  |
| Owner                 | Magentus Training Manager                        | Date Created      | 30/07/2022   |
| File Reference        | CRIS_CRIB_VR_410_Augnito_Voice_Commands_MAG_V2.0 |                   |  |
| Product Version       | 2.12.08 onwards                                  |                   |  |
| <b>Change History</b> |  |                   |  |
| Issue                 | Date   | Author            | Details of Change                                    |
| D1.0                  | 30/07/2022                                       | Emma Brown        | Reformat   |
| D1.1                  | 27/09/2022                                       | Emma Brown        | Corrections  |
| V1.0                  | 20/10/2022                                       | Emma Savage-Mady  | Updated in conjunction with validated commands list. |
| V1.1                  | 22/02/2023                                       | Emma Savage-Mady  | Updated to add Paragraph command.                    |
| V2.0                  | 10/11/2023                                       | Christine Anthony | Rebranded to Magentus.                               |
|                       |  |                   |  |
|                       |  |                   |  |
|                       |  |                   |  |
|                       |  |                   |  |
|                       |  |                   |  |
|                       |  |                   |  |
| Reviewed by           |  | Date              |  |
| Authorised by         |  | Date              |  |
| Review Date           |  |                   |  |