



core CRIS

Release 2.12.15

Release Notes (External)



Table of Contents

Copyright	3
Confidentiality	3
Document Control	3
1 Introduction	4
1.1 Purpose	4
1.2 Audience	4
1.3 Previous Release	4
2 Functionality added in this release.	5
2.1 RAD-34438 Medica Addendum Issues.	5
3 Fixes to existing functionality in this release	6
4 Known issues.	7
5 System Requirements and Supported Platforms	8



Copyright

© Magentus 2024

Registered Office: 3rd Floor i2 Mansfield, Hamilton Court, Oakham Business Park, Mansfield, NG18 5FB

These materials are or comprise restricted and proprietary confidential information of Magentus. Disclosure to or use by the recipient shall not convey any intellectual property rights in these materials. The right to use these materials by the recipient is subject to restrictions and limitations contained in the Provision of Integrated Care Record System and Associated Services Agreement and related agreements.

Confidentiality

All information in this document is provided in confidence for the sole purpose of adjudication of the document and shall not be used for any other purpose and shall not be published or disclosed wholly or in part to any other party without Magentus prior permission in writing and shall be held in safe custody. These obligations shall not apply to information which is published or becomes known legitimately from some source other than Magentus. Many of the product, service and company names referred to in this document are trademarks or registered trademarks. They are all hereby acknowledged.

Document Control

Title	core CRIS 2.12.15 External Release Notes		
Owner	CTO	Date Created	16/10/2024
File Reference	CRIS External Release Notes v2.12.15.docx		
Change History			
Issue	Date	Author	Details of Change
D1.0	16/10/2024	Gary Glover	First draft
v1.0	16/10/2024	Danielle West	Reviewed
Review Date			



1 Introduction

1.1 Purpose

This document will outline the functional changes made to this release. It will detail any new features, enhancements to existing functionality and list changes made due to existing functionality not working as expected.

1.2 Audience

This document is primarily for external stakeholders, so they have a clear understanding of what has changed within the application's functionality, compared to previous releases.

It is intended that the content of this document could be shared with external customers, but that most such communication would be through training materials and help files which are beyond the scope of this document.

1.3 Previous Release

The previous release of this product was as follows:

Version:	2.12.14
Date:	21/Feb/24



2 Functionality added in this release.

The following features and enhancements are included in this release:

2.1 RAD-34438 Medica Addendum Issues.

Jira Reference [CRIS-8346](#)

Currently, when we receive a verified report from an external source, the resulting verify triggers are suspended for a minimum of 2-minutes. During this delay, changes could be made to this report and the changes verified. Therefore, when the delay has elapsed, the status of the report is now at 'Verified Addendum (C)', the original verified report is overwritten. Some downstream PACS needs to capture the original report.

Resolution

Firstly, the '**REPORT.SendReportDelaySeconds**' XR setting **MUST** be set to **0** at the appropriate level. It defaults to 120 seconds if not explicitly set. So, set this at trust level.

This will ensure that all verified reports are sent downstream to PACS instantly. Please be aware, these 0 delays are for the verified reports only, there will still be a minimum of 120 seconds for unverified reports should they be configured to be sent to PACS.



3 Fixes to existing functionality in this release

There are no fixes in this release.



4 Known issues.

There were no new issues identified during testing of this release.



5 System Requirements and Supported Platforms

See attached Magentus Warranted Environment Specification v2.6



Magentus
Warranted Environn